


# LOCTEC<sup>®</sup>

Straight Plates 3.5 & 4.5  
Surgical Technique



Locking Compression Technology by aap

## Disclaimer

This surgical technique is exclusively intended for medical professionals, especially physicians, and therefore may not be regarded as a source of information for non-medical persons. The description of this surgical technique does not constitute medical advice or medical recommendations nor does it convey any diagnostic or therapeutic information on individual cases. Therefore, the attending physician is fully responsible for providing medical advice to the patient and obtaining the informed consent of the patient which this surgical technique does not supersede.

The description of this surgical technique has been compiled by medical experts and trained staff of aap Implantate AG with utmost diligence and to the best of their knowledge. However, aap Implantate AG excludes any liability for the completeness, accuracy, currentness, and quality of the information as well as for material or immaterial damages arising from the use of this information.

|   |           |
|---|-----------|
| <b>Introduction</b>                                     | <b>2</b>  |
| • Material  | 2         |
| • Indications/Contraindications                         | 3         |
| • Processing (Sterilization & Cleaning)                 | 3         |
| • MRI Safety Information                                | 3         |
| • Implants Small Fragment 3.5                           | 4         |
| • Implants Large Fragment 4.5                           | 4         |
| <b>Surgical Technique</b>                               | <b>6</b>  |
| • Preoperative planning                                 | 6         |
| • Preparing the plate                                   | 6         |
| • Reduction and primary fixation                        | 6         |
| • Insertion of cortical screws (gold)                   | 7         |
| • Insertion of cortical screws (blue)                   | 8         |
| • Insertion of locking screws (red) without compression | 9         |
| • Insertion of locking screws (red) with compression    | 10        |
| <b>Explantation</b>                                     | <b>12</b> |
| <b>Assembly instructions</b>                            | <b>13</b> |
| <b>Implants &amp; Instruments Small Fragment 3.5</b>    | <b>17</b> |
| <b>Implants &amp; Instruments Large Fragment 4.5</b>    | <b>19</b> |
| <b>Clinical Cases</b>                                   | <b>24</b> |

The LOQTEQ® Straight Plates 3.5 and 4.5 system comprises bone plates and screws as well as the necessary accessories for insertion. The flexibility of the system allows for safe reduction and stable fixation of various fracture patterns with respect to the indications mentioned in this technical guide. The design features of the implants in combination with standard surgical techniques and proven stability result in reliable constructs to support the bone during the healing process and allow for early mobilization of the patient.

## Material

The LOQTEQ® implants and instruments are manufactured using high-quality materials, which have been proven to be successful in medical technology for decades. The plates are made of pure titanium and the screws of a titanium alloy.

All materials employed comply with national and international standards. They are characterized by good biocompatibility, a high degree of safety against allergic reactions and good mechanical properties. LOQTEQ® implants show an excellent, highly polished surface.

## Indications/Contraindications

### Indications

#### LOQTEQ® Straight Plate 3.5

Fixation of fractures, osteotomies, and non-unions of the scapula, olecranon, humerus, radius, ulna, distal tibia, and fibula, particularly in osteoporotic bone.

#### LOQTEQ® Reconstruction Plate 3.5

Fixation of fractures, osteotomies, and non-unions of the scapula, olecranon, humerus, radius, ulna, distal tibia, and fibula, particularly in osteoporotic bone.

#### LOQTEQ® One Third Tubular Plate 3.5

Fixation of fractures, osteotomies, and non-unions of the scapula, olecranon, humerus, radius, ulna, distal tibia, and fibula, particularly in osteoporotic bone.

#### LOQTEQ® Narrow Plate 4.5 and LOQTEQ® Broad Plate 4.5

Fixation of various long bones, such as the humerus, femur and tibia. It is also for use in fixation of osteopenic bone and fixation of non-unions or malunions.

### Contraindications

- Infection or inflammation (local or systemic)
- Allergies to the implant material
- Acute and chronic osteomyelitis at or close to the surgical field
- High anesthesia risk patients
- Severe soft tissue swelling impacting normal wound healing
- Insufficient soft tissue coverage
- Fractures in children and adolescents with epiphyseal plates not yet ossified

#### ◆ Caution:

aap bone screws are neither designed nor approved for bolting or fixation of any elements (pediculi) of the cervical, thoracic or lumbar spine.

Detailed information on indications, contraindications and a complete list of adverse effects is included in the instructions for use.

## Processing (Sterilization & Cleaning)

aap markets unsterilized products which are appropriately labeled and must be appropriately processed before use (see Instructions for Use, chapter "Processing of Medical Devices").

Never use damaged implants or implants from damaged packaging.

## MRI Safety Information

The LOQTEQ® implant systems have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of LOQTEQ® implant systems in the MR environment is unknown.

Scanning a patient who has this medical device may result in patient injury.

## Implantats Small Fragment 3.5

Straight Plate 3.5, locking-compression hole



Reconstruction Plate 3.5, round locking hole



1/3 Tubular Plate 3.5, round locking hole



### LOQTEQ® Cortical Screw 3.5

- ▶ locking
- ▶ self-tapping



- ▶ red
- ▶ T15



Gliding-  
locking hole



For use in locking-compression hole:  
with/without compression.

- ▶ Drill (blue/red):  
ø2.7

### LOQTEQ® Cortical Screw 3.5, small Head

- ▶ locking
- ▶ self-tapping



- ▶ blue
- ▶ T15



Round-  
locking hole



For use in round locking  
hole.  
**Small head – do not use in  
gliding-locking hole!**

- ▶ Drill (blue/red):  
ø2.7

### Cortical Screw 3.5

- ▶ non-locking
- ▶ self-tapping



- ▶ gold
- ▶ T15



All plate holes

For use in all plate holes:  
with/without compression  
and as lag screw.

- ▶ Drill (blue/red):  
ø2.7



## Implantats Large Fragment 4.5

Straight plate 4.5 narrow, locking-compression hole



Straight plate 4.5 broad, locking-compression hole



### LOQTEQ® Cortical Screw 4.5

- ▶ locking
- ▶ self-tapping



- ▶ red
- ▶ T25



Gliding-  
locking hole



For use in locking-compression hole:  
with/without compression.

- ▶ Drill (blue/red):  
ø3.8



### Cortical Screw 4.5

- ▶ non-locking
- ▶ self-tapping




- ▶ gold
- ▶ T25



All plate holes

For use in all plate holes:  
with/without compression  
and as lag screw.

- ▶ Drill:  
ø3.2



The LOQTEQ® System includes many different plate types for use in the treatment of a wide number of fracture types. The primary focus of this surgical technique is the application of straight plates and does not deal with any specific fracture type. Please refer to special literature for more specific fracture treatment options.

The combination of angular stability and fracture compression in the novel gliding hole offers the following options of screw fixation:

- Locked fracture compression with locking screw
- Fracture compression with non-locking screw
- Locking and non-locking screws in neutral position

## Surgical Technique

### INSTRUMENTS SMALL FRAGMENT

Bending iron 1 for small fragment plates, closed  
Bending iron 2 for small fragment plates, closed  
Drill guide for gliding hole LOQTEQ® 3.5, I-ø 2.8, red  
Drill guide for round hole LOQTEQ® 3.5, I-ø 2.8, blue  
Reduction sleeve for K-wire ø1.6  
K-wire with trocar point, ø1.6, L 150

### ART.-NO.


IP 8405-00  
IP 8405-50  
IU 8166-10  
IU 8166-20  
IU 8166-16  
NK 0016-15

### INSTRUMENTS LARGE FRAGMENT

Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm  
Reduction sleeve for K-wire ø2.0  
K-wire with trocar point, ø2.0, L 250

### ART.-NO.

IU 8167-10  
IU 8167-15  
NK 0020-25



### Preoperative planning

- Evaluate the fracture situation and select the appropriate plate size and position with an X-ray. Consider the use of independent lag screws, if necessary.

### Preparing the Plate

- Select the plate that fits the fracture pattern and patient's anatomy.
- **NOTE:**  
Anatomically pre-contoured plates minimize the need for intraoperative bending. If necessary, the plates may be contoured.

#### ◆ CAUTION:

When contouring implants to the given shape of the bones, the implants should not be bent repeatedly as this may result in implant failure. Do not distort the implants or bend them within the confines of the plate holes and avoid any sharp-edged damage by instruments. Angle-stable plates are applied according to the concept of "internal fixation". The implants support bone healing through bone fixation but are not designed to withstand continuous loading without bone support and fracture healing. Since they are subject to material fatigue, medical devices applied in fracture fixation cannot be expected to withstand the same activities as normal healthy bone. Therefore, this support system for fracture healing is not as stable, strong and resilient as normal human bone.

### Reduction and primary fixation

- Reduce and temporarily secure the fracture. Care must be taken when positioning K-wires or independent lag screws, that they do not interfere with the later plate position.
- Insert and position the plate over the fracture site. The plate is fixed to the bone with K-wires. Confirm anatomic reduction and plate position using fluoroscopy.



## Insertion of cortical screws (gold)



### INSTRUMENTS SMALL FRAGMENT

Double drill guide  $\varnothing 2.7/3.5$ , with spring aided centering  
 Twist drill  $\varnothing 2.7$ , L 150, coil 50, quick coupling, single use  
 Twist drill  $\varnothing 2.7$ , L 150, coil 50, quick coupling, scaled, single use  
 Twist drill  $\varnothing 3.5$ , L 110, coil 50, quick coupling, single use  
 Depth gauge for screws  $\varnothing 3.5-4.0$ , up to L 90  
 Screwdriver Duo, T15, quick coupling  
 Large handle, cannulated, quick coupling

### ART.-NO.

IU 8116-60  
 IU 7427-15-1U  
 IU 7427-16-1U  
 IU 7435-00-1U  
 IS 7904-20  
 IU 7825-56  
 IU 7706-00

### INSTRUMENTS LARGE FRAGMENT

Double drill guide  $\varnothing 3.2/4.5$ , with spring aided centering  
 Twist drill  $\varnothing 3.2$ , L 195, coil 50, quick coupling, single use  
 Twist drill  $\varnothing 4.5$ , L 145, coil 50, quick coupling, single use  
 Depth gauge for screws  $\varnothing 4.5 - 6.5$ , up to L 100  
 Screwdriver Duo, T25, quick coupling  
 Large handle, cannulated, quick coupling

### ART.-NO.

IU 8117-50  
 IU 7432-30-1U  
 IU 7445-00-1U  
 IS 7905-20  
 IU 7835-56  
 IU 7706-00



### ◆ CAUTION:

Cortical screws are offered with different head diameters. Follow the instructions on page 15 for choosing the appropriate screw!

### • NOTE:

If a combination of locking and non-locking screws is used, non-locking cortical screws must be inserted first.

- To insert a cortical screw (gold), place the double drill guide in the plate hole and press it down. Choose an appropriate drill and drill through both cortices. Determine the length of the screw using the depth gauge and insert a screw of appropriate length using the screwdriver.

### • NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Check plate position using fluoroscopy and adjust if required.
- When using a cortical screw as a lag screw, start with the drill for the gliding hole and the appropriate side of the double drill guide. Drill through the near cortex or perforating the fracture line, center the other side of the drill guide in the gliding hole and drill the core hole with an appropriate diameter drill through the far cortex. Determine the screw length using the depth gauge and insert a non-locking coral screw of the appropriate length.



## only for Reconstruction Plate 3.5 and 1/3 Tubular Plate 3.5

Insertion of cortical screws  
(blue)

## INSTRUMENTS

Drill guide for round hole LOQTEQ® 3.5, I-ø 2.8, blue  
 Twist drill ø2.7, L 150, coil 50, quick coupling, single use  
 Twist drill ø2.7, L 150, coil 50, quick coupling, scaled, single use  
 Stop ring for depth measurement, SF  
 Depth gauge for screws ø3.5–4.0, up to L 90  
 Screwdriver Duo, T15, quick coupling  
 Handle with quick coupling, with torque limiter, 2.0Nm

## ART.-NO.

IU 8166-20  
 IU 7427-15-1U  
 IU 7427-16-1U  
 IU 8166-06  
 IS 7904-20  
 IU 7825-56  
 IU 7707-20



3.5



- Insert a drill guide ø2.7 mm (blue) into any chosen round hole and drill to the desired depth using a drill (blue/red).

## ◆ CAUTION:

The screwdriver duo is not intended for screwing the drill guide into the plate.

- The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.
- The stop ring facilitates reading off the calibration when attached to the drill. Push it down to the guide sleeve and remove the drill for reading the drilling depth in the gap of the ring.

## • NOTE:

The screwdriver duo facilitates manual removal of the drill guide.

- Select a locking screw (blue) of the proper length. Loosely insert the screw using the screwdriver T15 manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface.

## • NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Finish the screw manually using the screwdriver T15 with the torque limiting handle 2.0 Nm. Optimal locking should be achieved with an audible and tactile click of the torque limiter.

## ◆ CAUTION:

As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.

- Follow these instructions to insert further screws in the plate holes depending on the fracture pattern. Finally, confirm that all screw heads are flush with the plate surface. Check the result using fluoroscopy and adjust screw positioning or length as necessary.

# Insertion of locking compression screws (red) without compression



## INSTRUMENTS SMALL FRAGMENT

Drill guide for gliding hole LOQTEQ® 3.5, I-ø 2.8, red  
Twist drill ø2.7, L 150, coil 50, quick coupling, single use  
Twist drill ø2.7, L 150, coil 50, quick coupling, scaled, single use  
Stop ring for depth measurement, SF  
Depth gauge for screws ø3.5-4.0, up to L 90  
Screwdriver Duo, T15, quick coupling  
Handle with quick coupling, with torque limiter, 2.0Nm

## ART.-NO.

IU 8166-10  
IU 7427-15-1U  
IU 7427-16-1U  
IU 8166-06  
IS 7904-20  
IU 7825-56  
IU 7707-20

## INSTRUMENTS LARGE FRAGMENT

Drill guide for gliding hole LOQTEQ® 4.5, I-ø 3.9, red  
Twist drill ø3.8, L 180, coil 50, quick coupling, single use  
Stop ring for depth measurement, LF  
Depth gauge for screws ø4.5-6.5, up to L 100  
Screwdriver Duo, T25, quick coupling  
Handle with quick coupling, with torque limiter 3.5Nm

## ART.-NO.

IU 8167-10  
IU 7438-18-1U  
IU 8184-03  
IS 7905-20  
IU 7835-56  
IU 7707-35



- Insert a drill guide (red) into any chosen gliding hole and drill to the desired depth using a drill (blue/red).

### ◆ CAUTION:

The screwdriver duo is not intended for screwing the drill guide into the plate.

- The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.
- The stop ring facilitates reading off the calibration when attached to the drill. Push it down to the guide sleeve and remove the drill for reading the drilling depth in the gap of the ring.

### • NOTE:

The screwdriver duo facilitates manual removal of the drill guide.

- Select a locking screw (rot) of the proper length. Loosely insert the screw using the screwdriver manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface.

### • NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Finish the screw manually using the screwdriver bit with the torque limiting handle. Optimal locking should be achieved with an audible and tactile click of the torque limiter.

### ◆ CAUTION:

As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.

### Insertion of locking compression screws (red) with compression



#### INSTRUMENTS SMALL FRAGMENT

Basic Insert for Load Drill Guide LOQTEQ® 3.5  
Load Drill guide LOQTEQ® 3.5, compression 1mm  
Load Drill guide LOQTEQ® 3.5, compression 2mm  
Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm

#### ART.-NO.

IU 8166-05  
IU 8166-01  
IU 8166-02  
IU 8166-03

#### INSTRUMENTS LARGE FRAGMENT

Basic insert for load drill guide LOQTEQ® 4.5  
Load Drill guide LOQTEQ® 4.5, compression 1mm  
Load Drill guide LOQTEQ® 4.5, compression 2mm  
Load drill guide LOQTEQ® 4.5, adjustable up to 2mm

#### ART.-NO.

IU 8167-05  
IU 8167-01  
IU 8167-02  
IU 8167-03





- LOQTEQ® gliding holes allow for fracture compression with subsequent locking fixation in one step. Load drill guides enable compressing fracture gaps of up to 2 mm.
- Screw the basic insert for load drill guides into a shaft hole near the fracture line or, if necessary, above the fracture line. Choose a load drill guide in accordance with the compression distance (1 mm or 2 mm), slide it on the basic insert and place the drill guide in the next plate hole, away from the fracture gap. Avoid pressure on the drill guide.

#### ◆ CAUTION:

The screwdriver duo is not intended for screwing the basic insert into the plate.

- If available, use the adjustable load drill guide. The fracture gap serves as orientation in setting the compression distance (max. 2 mm) by turning the wheel of the load drill guide until an appropriate gap opens in the upper part of the instrument.
- **NOTE:**  
Care should be taken when selecting the proper compression distance (1 mm or 2 mm). Avoid overcompression to ensure full locking of the screw, especially in hard bone.





#### INSTRUMENTS SMALL FRAGMENT

Twist drill ø2.7, L 150, coil 50, quick coupling, single use  
Twist drill ø2.7, L 150, coil 50, quick coupling, scaled, single use  
Depth gauge for screws ø3.5-4.0, up to L 90  
Screwdriver Duo, T15, quick coupling  
Handle with quick coupling, with torque limiter, 2.0 Nm  
Large handle, cannulated, quick coupling

#### ART.-NO.

IU 7427-15-1U  
IU 7427-16-1U  
IS 7904-20  
IU 7825-56  
IU 7707-20  
IU 7706-00

#### INSTRUMENTS LARGE FRAGMENT



Twist drill ø3.8, L 250, coil 50, quick coupling, single use  
Depth gauge for screws ø4.5-6.5, up to L 100  
Screwdriver Duo, T25, quick coupling  
Handle with quick coupling, with torque limiter 3.5 Nm  
Large handle, cannulated, quick coupling

#### ART.-NO.

IU 7438-25-1U  
IS 7905-20  
IU 7835-56  
IU 7707-35  
IU 7706-00



- Drill to the desired depth using a drill (blue/red) and remove the basic insert. The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.
- **NOTE:**  
The screwdriver duo facilitates manual removal of the basic insert.
- Select a locking screw (red) of the proper length. Loosely insert the screw using the screwdriver manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface. Finish the screw manually using the screwdriver bit with the torque limiting handle. Optimal locking should be achieved with an audible and tactile click of the torque limiter.
- ◆ **CAUTION:**  
As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.
- Alternatively, insert a non-locking cortical screw (gold) as a compression screw by placing the double drill guide without pressure in an off-center position in the plate hole. For insertion of the screw, follow the instructions on page 7.
- Follow these instructions to insert further screws in the plate holes depending on the fracture pattern. Finally, confirm that all screw heads are flush with the plate surface. Check the result using fluoroscopy and adjust screw positioning or length as necessary.



**INSTRUMENTS**


Explantation screwdriver T15, round handle  
Explantation screwdriver T25, round handle

**ART.-NO.**

IU 7811-15  
IU 7811-25

- **NOTE:**  
The screwdrivers T15 (IU 7825-56) and T25 (IU 7835-56) in the set are self-retaining and should not be used for screw explantation.
- Use the corresponding explantation screwdriver for safe removal of a screw. Explantation screwdrivers are not self-retaining, penetrate further into the screw head and thus permit a higher torque when removing screws. They are not included in the set as standard and must be ordered separately.
- Place an incision on the old scar. Manually undo all screws and sequentially remove them. After manually unlocking all screws, removal may be performed using a power tool.


## Disassembly



### Adjustable load drill guide (IU 8166-03 / IU 8167-03)

- Remove screws (item 4.1 and 4.2) using a hexagonal screwdriver 2.5
- Unscrew the set screw (item 3)
- Pull the compression block apart (items 1 and 2)

## Assembly



- Fit together the compression block (items 1 and 2)
- Insert the set screw (item 3) into the compression block, middle hole
- Insert the retaining screws (items 4.1 and 4.2) using a hexagonal screwdriver 2.5





### LOQTEQ® Straight Plate 3.5

| HOLES | LENGTH<br>(mm) | ART.-NO.     |
|-------|----------------|--------------|
| 4     | 60             | PG 3555-04-2 |
| 5     | 73             | PG 3555-05-2 |
| 6     | 86             | PG 3555-06-2 |
| 7     | 99             | PG 3555-07-2 |
| 8     | 112            | PG 3555-08-2 |
| 9     | 125            | PG 3555-09-2 |
| 10    | 138            | PG 3555-10-2 |
| 11    | 151            | PG 3555-11-2 |
| 12    | 164            | PG 3555-12-2 |
| 14    | 190            | PG 3555-14-2 |



### LOQTEQ® Reconstruction Plate 3.5

| HOLES | LENGTH<br>(mm) | ART.-NO.     |
|-------|----------------|--------------|
| 5     | 58             | PR 3550-05-2 |
| 6     | 70             | PR 3550-06-2 |
| 7     | 82             | PR 3550-07-2 |
| 8     | 94             | PR 3550-08-2 |
| 9     | 106            | PR 3550-09-2 |
| 10    | 118            | PR 3550-10-2 |
| 12    | 142            | PR 3550-12-2 |
| 14    | 166            | PR 3550-14-2 |
| 16    | 190            | PR 3550-16-2 |
| 18    | 214            | PR 3550-18-2 |



### LOQTEQ® 1/3 Tubular Plate 3.5

| HOLES | LENGTH<br>(mm) | ART.-NO.     |
|-------|----------------|--------------|
| 3     | 40             | PG 3553-03-2 |
| 4     | 52             | PG 3553-04-2 |
| 5     | 64             | PG 3553-05-2 |
| 6     | 76             | PG 3553-06-2 |
| 7     | 88             | PG 3553-07-2 |
| 8     | 100            | PG 3553-08-2 |
| 9     | 112            | PG 3553-09-2 |
| 10    | 124            | PG 3553-10-2 |
| 11    | 136            | PG 3553-11-2 |
| 12    | 148            | PG 3553-12-2 |



**LOQTEQ® Cortical Screw 3.5, T15,  
self-tapping**


|      |              |
|------|--------------|
| L 12 | SK 3525-12-2 |
| L 14 | SK 3525-14-2 |
| L 16 | SK 3525-16-2 |
| L 18 | SK 3525-18-2 |
| L 20 | SK 3525-20-2 |
| L 22 | SK 3525-22-2 |
| L 24 | SK 3525-24-2 |
| L 26 | SK 3525-26-2 |
| L 28 | SK 3525-28-2 |
| L 30 | SK 3525-30-2 |
| L 32 | SK 3525-32-2 |
| L 34 | SK 3525-34-2 |
| L 36 | SK 3525-36-2 |
| L 38 | SK 3525-38-2 |
| L 40 | SK 3525-40-2 |
| L 42 | SK 3525-42-2 |
| L 44 | SK 3525-44-2 |
| L 45 | SK 3525-45-2 |
| L 50 | SK 3525-50-2 |
| L 55 | SK 3525-55-2 |
| L 60 | SK 3525-60-2 |
| L 65 | SK 3525-65-2 |
| L 70 | SK 3525-70-2 |

**LOQTEQ® Cortical Screw 3.5,  
small head, T15, self-tapping**


|      |              |
|------|--------------|
| L 10 | SK 3526-10-2 |
| L 12 | SK 3526-12-2 |
| L 14 | SK 3526-14-2 |
| L 16 | SK 3526-16-2 |
| L 18 | SK 3526-18-2 |
| L 20 | SK 3526-20-2 |
| L 22 | SK 3526-22-2 |
| L 24 | SK 3526-24-2 |
| L 26 | SK 3526-26-2 |
| L 28 | SK 3526-28-2 |
| L 30 | SK 3526-30-2 |
| L 32 | SK 3526-32-2 |
| L 34 | SK 3526-34-2 |
| L 36 | SK 3526-36-2 |
| L 38 | SK 3526-38-2 |
| L 40 | SK 3526-40-2 |
| L 42 | SK 3526-42-2 |
| L 45 | SK 3526-45-2 |
| L 50 | SK 3526-50-2 |
| L 55 | SK 3526-55-2 |
| L 60 | SK 3526-60-2 |
| L 65 | SK 3526-65-2 |
| L 70 | SK 3526-70-2 |
| L 75 | SK 3526-75-2 |
| L 80 | SK 3526-80-2 |
| L 85 | SK 3526-85-2 |
| L 90 | SK 3526-90-2 |

**Cortical Screw 3.5,  
T15, self-tapping**


|      |              |
|------|--------------|
| L 10 | SK 3514-10-2 |
| L 12 | SK 3514-12-2 |
| L 14 | SK 3514-14-2 |
| L 16 | SK 3514-16-2 |
| L 18 | SK 3514-18-2 |
| L 20 | SK 3514-20-2 |
| L 22 | SK 3514-22-2 |
| L 24 | SK 3514-24-2 |
| L 26 | SK 3514-26-2 |
| L 28 | SK 3514-28-2 |
| L 30 | SK 3514-30-2 |
| L 32 | SK 3514-32-2 |
| L 34 | SK 3514-34-2 |
| L 36 | SK 3514-36-2 |
| L 38 | SK 3514-38-2 |
| L 40 | SK 3514-40-2 |
| L 42 | SK 3514-42-2 |
| L 45 | SK 3514-45-2 |
| L 50 | SK 3514-50-2 |
| L 55 | SK 3514-55-2 |
| L 60 | SK 3514-60-2 |
| L 65 | SK 3514-65-2 |
| L 70 | SK 3514-70-2 |
| L 75 | SK 3514-75-2 |
| L 80 | SK 3514-80-2 |
| L 85 | SK 3514-85-2 |
| L 90 | SK 3514-90-2 |

### Cancellous Screw 4.0, small head, T15



|      |       |              |
|------|-------|--------------|
| L 10 | TL 5  | SP 4030-10-2 |
| L 12 | TL 5  | SP 4030-12-2 |
| L 14 | TL 5  | SP 4030-14-2 |
| L 16 | TL 6  | SP 4030-16-2 |
| L 18 | TL 7  | SP 4030-18-2 |
| L 20 | TL 8  | SP 4030-20-2 |
| L 22 | TL 9  | SP 4030-22-2 |
| L 24 | TL 10 | SP 4030-24-2 |
| L 26 | TL 12 | SP 4030-26-2 |
| L 28 | TL 14 | SP 4030-28-2 |
| L 30 | TL 14 | SP 4030-30-2 |
| L 32 | TL 14 | SP 4030-32-2 |
| L 34 | TL 14 | SP 4030-34-2 |
| L 36 | TL 14 | SP 4030-36-2 |
| L 38 | TL 14 | SP 4030-38-2 |
| L 40 | TL 14 | SP 4030-40-2 |
| L 42 | TL 15 | SP 4030-42-2 |
| L 45 | TL 15 | SP 4030-45-2 |
| L 50 | TL 15 | SP 4030-50-2 |
| L 55 | TL 16 | SP 4030-55-2 |
| L 60 | TL 16 | SP 4030-60-2 |
| L 65 | TL 16 | SP 4030-65-2 |
| L 70 | TL 16 | SP 4030-70-2 |
| L 75 | TL 16 | SP 4030-75-2 |
| L 80 | TL 16 | SP 4030-80-2 |
| L 85 | TL 16 | SP 4030-85-2 |
| L 90 | TL 16 | SP 4030-90-2 |

### Cancellous Screw 4.0, kl. Kopf, T15, full thread



|      |              |
|------|--------------|
| L 10 | SP 4035-10-2 |
| L 12 | SP 4035-12-2 |
| L 14 | SP 4035-14-2 |
| L 16 | SP 4035-16-2 |
| L 18 | SP 4035-18-2 |
| L 20 | SP 4035-20-2 |
| L 22 | SP 4035-22-2 |
| L 24 | SP 4035-24-2 |
| L 26 | SP 4035-26-2 |
| L 28 | SP 4035-28-2 |
| L 30 | SP 4035-30-2 |
| L 32 | SP 4035-32-2 |
| L 34 | SP 4035-34-2 |
| L 36 | SP 4035-36-2 |
| L 38 | SP 4035-38-2 |
| L 40 | SP 4035-40-2 |
| L 42 | SP 4035-42-2 |
| L 45 | SP 4035-45-2 |
| L 50 | SP 4035-50-2 |
| L 55 | SP 4035-55-2 |
| L 60 | SP 4035-60-2 |
| L 65 | SP 4035-65-2 |
| L 70 | SP 4035-70-2 |
| L 75 | SP 4035-75-2 |
| L 80 | SP 4035-80-2 |
| L 85 | SP 4035-85-2 |

### Washer

I-ø 4.4 mm, A-ø 8.0 mm, Titanium



SU 0448-00-2




Bending iron 1 for small fragment plates, closed  
Bending iron 2 for small fragment plates, closed

IP 8405-00  
IP 8405-50



Depth gauge for screws ø3.5 - 4.0, up to L 90

IS 7904-20



Twist drill ø2.7, L 150, coil 50, quick coupling, single use  
Twist drill ø2.7, L 150, coil 50, quick coupling, scaled, single use  
Twist drill ø2.7, L 220, coil 50, quick coupling, single use  
Twist drill ø2.7, L 220, coil 50, quick coupling, scaled, single use  
Twist drill ø3.5, L 110, coil 50, quick coupling, single use

IU 7427-15-1U  
IU 7427-16-1U  
IU 7427-22-1U  
IU 7427-23-1U  
IU 7435-00-1U



Large handle, cannulated, quick coupling

IU 7706-00



Handle with quick coupling, with torque limiter 2.0Nm

IU 7707-20



Screwdriver Duo, T15, quick coupling

IU 7825-56 ★



Double drill guide  $\varnothing 2.7/3.5$ , with spring aided centering

IU 8116-60 ★



Load Drill guide LOQTEQ® 3.5, compression 1mm

IU 8166-01

Load Drill guide LOQTEQ® 3.5, compression 2mm

IU 8166-02

Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm

IU 8166-03

Basic insert for load drill guide LOQTEQ® 3.5

IU 8166-05



Stop ring for depth measurement, SF

IU 8166-06



Drill guide for gliding hole LOQTEQ® 3.5, I- $\varnothing$  2.8, red

IU 8166-10




Reduction sleeve for K-wire  $\varnothing 1.6$

IU 8166-16



Drill guide for round hole LOQTEQ® 3.5, I- $\varnothing$  2.8, blue

IU 8166-20



Caddy for K-wire L 150

IC 0006-15

K-wire with trocar point,  $\varnothing 1.6$ , L 150

NK 0016-15





### LOQTEQ® Narrow Plate 4.5

| HOLES | LENGTH<br>(mm) | ART.-NO.     |
|-------|----------------|--------------|
| 4     | 72             | PG 4555-04-2 |
| 5     | 90             | PG 4555-05-2 |
| 6     | 108            | PG 4555-06-2 |
| 7     | 126            | PG 4555-07-2 |
| 8     | 144            | PG 4555-08-2 |
| 9     | 162            | PG 4555-09-2 |
| 10    | 180            | PG 4555-10-2 |
| 11    | 198            | PG 4555-11-2 |
| 12    | 216            | PG 4555-12-2 |
| 14    | 252            | PG 4555-14-2 |
| 16    | 288            | PG 4555-16-2 |



### LOQTEQ® Broad Plate 4.5

| HOLES | LENGTH<br>(mm) | ART.-NO.     |
|-------|----------------|--------------|
| 6     | 115            | PG 4556-06-2 |
| 7     | 133            | PG 4556-07-2 |
| 8     | 150            | PG 4556-08-2 |
| 9     | 168            | PG 4556-09-2 |
| 10    | 186            | PG 4556-10-2 |
| 11    | 204            | PG 4556-11-2 |
| 12    | 222            | PG 4556-12-2 |
| 14    | 257            | PG 4556-14-2 |
| 16    | 293            | PG 4556-16-2 |
| 18    | 328            | PG 4556-18-2 |

**LOQTEQ® Cortical Screw 4.5,  
T25, self-tapping**



|      |              |
|------|--------------|
| L 14 | SK 4525-14-2 |
| L 16 | SK 4525-16-2 |
| L 18 | SK 4525-18-2 |
| L 20 | SK 4525-20-2 |
| L 22 | SK 4525-22-2 |
| L 24 | SK 4525-24-2 |
| L 26 | SK 4525-26-2 |
| L 28 | SK 4525-28-2 |
| L 30 | SK 4525-30-2 |
| L 32 | SK 4525-32-2 |
| L 34 | SK 4525-34-2 |
| L 36 | SK 4525-36-2 |
| L 38 | SK 4525-38-2 |
| L 40 | SK 4525-40-2 |
| L 42 | SK 4525-42-2 |
| L 45 | SK 4525-45-2 |
| L 50 | SK 4525-50-2 |
| L 55 | SK 4525-55-2 |
| L 60 | SK 4525-60-2 |
| L 65 | SK 4525-65-2 |
| L 70 | SK 4525-70-2 |
| L 75 | SK 4525-75-2 |
| L 80 | SK 4525-80-2 |
| L 85 | SK 4525-85-2 |
| L 90 | SK 4525-90-2 |

**Cortical Screw 4.5,  
T25, self-tapping**



|      |              |
|------|--------------|
| L 20 | SK 4514-20-2 |
| L 22 | SK 4514-22-2 |
| L 24 | SK 4514-24-2 |
| L 26 | SK 4514-26-2 |
| L 28 | SK 4514-28-2 |
| L 30 | SK 4514-30-2 |
| L 32 | SK 4514-32-2 |
| L 34 | SK 4514-34-2 |
| L 36 | SK 4514-36-2 |
| L 38 | SK 4514-38-2 |
| L 40 | SK 4514-40-2 |
| L 42 | SK 4514-42-2 |
| L 45 | SK 4514-45-2 |
| L 50 | SK 4514-50-2 |
| L 55 | SK 4514-55-2 |
| L 60 | SK 4514-60-2 |
| L 65 | SK 4514-65-2 |
| L 70 | SK 4514-70-2 |
| L 75 | SK 4514-75-2 |
| L 80 | SK 4514-80-2 |
| L 85 | SK 4514-85-2 |
| L 90 | SK 4514-90-2 |



Depth gauge for screws  $\varnothing 4.5$  - 6.5, up to L 100

IS 7905-20



Twist drill  $\varnothing 3.2$ , L 195, coil 50, quick coupling, single use  
Twist drill  $\varnothing 3.8$ , L 180, coil 50, quick coupling, single use  
Twist drill  $\varnothing 3.8$ , L 250, coil 50, quick coupling, single use  
Twist drill  $\varnothing 4.5$ , L 145, coil 50, quick coupling, single use

IU 7432-30-1U  
IU 7438-18-1U  
IU 7438-25-1U  
IU 7445-00-1U



Large handle, cannulated, quick coupling

IU 7706-00



Handle with quick coupling, with torque limiter 3.5Nm

IU 7707-35



Screwdriver Duo, T25, quick coupling

IU 7835-56 ★



Double drill guide  $\varnothing 3.2/4.5$ , with spring aided centering

IU 8117-50





Load drill guide LOQTEQ® 4.5, compression 1mm  
Load drill guide LOQTEQ® 4.5, compression 2mm  
Load drill guide LOQTEQ® 4.5, adjustable up to 2mm  
Basic insert for load drill guide LOQTEQ® 4.5

IU 8167-01  
IU 8167-02  
IU 8167-03  
IU 8167-05



Drill guide for gliding hole LOQTEQ® 4.5, I-ø 3.9, red

IU 8167-10




Reduction sleeve for K-wire ø2.0

IU 8167-15



Stop ring for depth measurement, LF

IU 8184-03




Caddy for K-wire L 250  
K-wire with trocar point, ø2.0, L 250

IC 0006-25  
NK 0020-25

Preoperative




Postoperative




Clinical case and CT images with the kind permission of  
Asklepios Clinic Weißenfels, Germany

## Ulna shortening osteotomy

Preoperative



Postoperative



Clinical case and  
CT images with the kind  
permission of  
University Hospital Gießen  
and Marburg, Germany

Subject to technical modifications,  
errors and misprints.

© aap Implantate AG  
WP 40P001 US / 2111-1

**aap Implants Inc.**  
260 Peachtree Street NW • Suite 2200  
Atlanta • GA 30303 USA

Phone +1 678-942-3791  
e-Fax +1 877-373-0637

[aapcs@gmlx.us.com](mailto:aapcs@gmlx.us.com)  
[www.aap-implants.com](http://www.aap-implants.com)

# LOQTEC®

---



aap Implants Inc.  
260 Peachtree Street NW • Suite 2200  
Atlanta • GA 30303 USA

Phone +1 678-942-3791  
e-Fax +1 877-373-0637

[aapcs@gmlx.us.com](mailto:aapcs@gmlx.us.com)  
[www.aap-implants.com](http://www.aap-implants.com)