


LOCTEC[®]

Straight Plates 3.5 & 4.5
Surgical Technique



Disclaimer

This surgical technique is exclusively intended for medical professionals, especially physicians, and therefore may not be regarded as a source of information for non-medical persons. The description of this surgical technique does not constitute medical advice or medical recommendations nor does it convey any diagnostic or therapeutic information on individual cases. Therefore, the attending physician is fully responsible for providing medical advice to the patient and obtaining the informed consent of the patient which this surgical technique does not supersede.

The description of this surgical technique has been compiled by medical experts and trained staff of aap Implantate AG with utmost diligence and to the best of their knowledge. However, aap Implantate AG excludes any liability for the completeness, accuracy, currentness, and quality of the information as well as for material or immaterial damages arising from the use of this information.

Introduction	2
• Material	2
• Intended use	2
• Indications/Contraindications	3
• Processing (Sterilization & Cleaning)	3
• MRI Safety Information	3
• Implants Small Fragment 3.5	4
• Implants Large Fragment 4.5	4
Surgical Technique	6
• Preoperative planning	6
• Preparing the plate	6
• Reduction and primary fixation	6
• Insertion of cortical screws (gold)	7
• Insertion of cortical screws (blue)	8
• Insertion of locking screws (red) without compression	9
• Insertion of locking screws (red) with compression	10
Explantation	12
Assembly instructions	13
Implants & Instruments Small Fragment 3.5	17
Implants & Instruments Large Fragment 4.5	19
Clinical Cases	24

The LOQTEQ® Straight Plates 3.5 and 4.5 system comprises bone plates and screws as well as the necessary accessories for insertion. The flexibility of the system allows for safe reduction and stable fixation of various fracture patterns with respect to the indications mentioned in this technical guide. The design features of the implants in combination with standard surgical techniques and proven stability result in reliable constructs to support the bone during the healing process and allow for early mobilization of the patient.

Material

The LOQTEQ® implants and instruments are manufactured using high-quality materials, which have been proven to be successful in medical technology for decades. The plates are made of pure titanium and the screws of a titanium alloy.

All materials employed comply with national and international standards. They are characterized by good biocompatibility, a high degree of safety against allergic reactions and good mechanical properties. LOQTEQ® implants show an excellent, highly polished surface.

Intended Use

LOQTEQ® Small Fragment Set 3.5

The plate and screw implants in the LOQTEQ® Small Fragment Set 3.5 are intended for the temporary fixation, correction or stabilization of the radius, ulna, humerus, scapula, fibula, pelvis and/or calcaneus. The implants are intended for single use in human bone.

LOQTEQ® Large Fragment Set 4.5

The plate and screw implants of the LOQTEQ® Large Fragment Set 4.5 are intended for the temporary fixation, correction or stabilization of diaphyseal fractures of the humerus, femur and/or tibia. The implants are intended for single use in human bone.

Indikationen/Kontraindikationen

Indications for Use

LOQTEQ® Straight Plate 3.5

Diaphyseal fractures of the radius; diaphyseal fractures of the ulna; diaphyseal fractures of the humerus; pseudarthrosis of the humeral diaphysis.

LOQTEQ® Reconstruction Plate 3.5

Displaced fractures of the scapula; distal fractures of the humerus AO type 13-A2, -A3, -B, and -C; olecranon fractures AO type 21-A1.3, -B1, -B3.1, -B3.3, and -C; pelvic fractures, in particular at the pubic symphysis and the inner pelvic girdle; acetabular fractures.

LOQTEQ® One Third Tubular Plate 3.5

Olecranon fractures; Weber fibular fractures; calcaneal fractures; fractures of the 1st metatarsal; metatarsal-phalangeal arthrodesis; LisFranc arthrodesis.

LOQTEQ® Narrow Plate 4.5 and LOQTEQ® Broad Plate 4.5

Diaphyseal fractures of the humerus; diaphyseal fractures of the femur; diaphyseal fractures of the tibia; fixation of periprosthetic fractures, non-unions and mal-unions.

Absolute Contraindications

- Infection or inflammation (local or systemic)
- Allergies to the implant material
- Acute or chronic osteomyelitis at or close to the surgical field
- Unacceptably high anesthesia risk
- Severe soft tissue swelling compromising normal wound healing
- Insufficient soft tissue coverage
- Fractures in children and adolescents with epiphyseal plates that are not yet ossified

Caution:

aap products are not approved for the spine.

Detailed information on indications, contraindications and a complete list of adverse effects is included in the instructions for use.

Processing (Sterilization & Cleaning)

aap markets unsterilized products which are appropriately labeled and must be appropriately processed before use (see Instructions for Use, chapter "Processing of products").

Never use damaged implants or implants from damaged packaging.

MR Safety Information

Non-clinical testing has demonstrated that the LOQTEQ® Small Fragment Set 3.5 and LOQTEQ® Large Fragment Set 4.5 systems are **MR Conditional**. Further information is included in the Instructions for Use that are enclosed with the products.



Implantats Small Fragment 3.5

Straight Plate 3.5, locking-compression hole



Reconstruction Plate 3.5, round locking hole



1/3 Tubular Plate 3.5, round locking hole



LOQTEQ® Cortical Screw 3.5

- ▶ locking
- ▶ self-tapping



- ▶ red
- ▶ T15




Gliding-
locking hole

LOQTEQ® Cortical Screw 3.5, small Head

- ▶ locking
- ▶ self-tapping



- ▶ blue
- ▶ T15



Round-
locking hole

Cortical Screw 3.5

- ▶ non-locking
- ▶ self-tapping



- ▶ gold
- ▶ T15




All plate holes

Cortical Screw 3.5

- ▶ non-locking
- ▶ self-tapping



- ▶ gold
- ▶ SW 2.5



large head for
gliding locking hole*

small head for
round locking hole*

For use in locking-compression hole: with/without compression.

- ▶ Drill (blue/red):
ø2.7

For use in round locking hole.
Small head – do not use in gliding-locking hole!

- ▶ Drill (blue/red):
ø2.7

For use in all plate holes: with/without compression and as lag screw.

- ▶ Drill (blue/red):
ø2.7

For use in all plate holes: with/without compression and as lag screw.

- ▶ Drill
ø2.5



* see Cautions on page 15

Implantats Large Fragment 4.5

Straight plate 4.5 narrow, locking-compression hole



Straight plate 4.5 broad, locking-compression hole



LOQTEQ® Cortical Screw 4.5

- ▶ locking
- ▶ self-tapping



- ▶ red
- ▶ T25



Gliding-
locking hole



For use in locking-compression hole:
with/without compression.

- ▶ Drill (blue/red):
ø3.8



Cortical Screw 4.5

- ▶ non-locking
- ▶ self-tapping



- ▶ gold
- ▶ T25



All plate holes

For use in all plate holes:
with/without compression
and as lag screw.

- ▶ Drill:
ø3.2




Cortical Screw 4.5

- ▶ non-locking
- ▶ self-tapping



- ▶ gold
- ▶ SW 3.5



For use in all plate holes:
with/without compression
and as lag screw.

- ▶ Drill:
ø3.2

The LOQTEQ® System includes many different plate types for use in the treatment of a wide number of fracture types. The primary focus of this surgical technique is the application of straight plates and does not deal with any specific fracture type. Please refer to special literature for more specific fracture treatment options.

The combination of angular stability and fracture compression in the novel gliding hole offers the following options of screw fixation:

- Locked fracture compression with locking screw
- Fracture compression with non-locking screw
- Locking and non-locking screws in neutral position

Surgical Technique

INSTRUMENTS SMALL FRAGMENT

Bending iron 1 for small fragment plates, closed
Bending iron 2 for small fragment plates, closed
Drill guide for gliding hole LOQTEQ® 3.5, I-ø 2.8, red
Drill guide for round hole LOQTEQ® 3.5, I-ø 2.8, blue
Reduction sleeve for K-wire ø1.6
K-wire with trocar point, ø1.6, L 150

ART.-NO.


IP 8405-00
IP 8405-50
IU 8166-10
IU 8166-20
IU 8166-16
NK 0016-15

INSTRUMENTS LARGE FRAGMENT

Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm
Reduction sleeve for K-wire ø2.0
K-wire with trocar point, ø2.0, L 250

ART.-NO.

IU 8167-10
IU 8167-15
NK 0020-25



Präoperative Planung

- Evaluate the fracture situation and select the appropriate plate size and position with an X-ray. Consider the use of independent lag screws, if necessary.

Preparing the Plate

- Select the plate that fits the fracture pattern and patient's anatomy.

◆ HINWEIS:

Anatomically pre-contoured plates minimize the need for intra-operative bending. If necessary, the plates may be contoured.

◆ CAUTION:

Anatomically preformed plates should not be bent where possible. If plates are adapted to anatomical bone structures, the implants should not be bent back and forth repeatedly and excessively as this may result in implant failure. Damage caused by sharp edges should be avoided when bending. Locking plates should in principle be bent in the area between the holes only. Bending plates along locking holes may impair or even abolish their function completely. If angular stability is compromised by bending, a non-locking screw should be used.

Reduction and primary fixation

- Reduce and temporarily secure the fracture. Care must be taken when positioning K-wires or independent lag screws, that they do not interfere with the later plate position.
- Insert and position the plate over the fracture site. The plate is fixed to the bone with K-wires. Confirm anatomic reduction and plate position using fluoroscopy.

Insertion of cortical screws (gold)



INSTRUMENTS SMALL FRAGMENT

	◆	★
Double drill guide with spring aided centering	IU 8116-50	IU 8116-60
Twist drill, quick coupling	IU 7425-00	IU 7427-15
Twist drill ø2.7, L 150, coil 50, quick coupling, scaled	–	IU 7427-16
Twist drill ø3.5, L 110, coil 50, quick coupling	IU 7435-00	IU 7435-00
Depth gauge for screws ø3.5-4.0, up to L 90	IS 7904-20	IS 7904-20
Screwdriver, quick coupling	IU 7825-00	IU 7825-56
Handle for quick coupling, large, cannulated	IU 7706-00	IU 7706-00

INSTRUMENTS LARGE FRAGMENT

	◆	★
Double drill guide ø3.2/4.5, with spring aided centering	IU 8117-50	IU 8117-50
Twist drill ø3.2, L 195, coil 50, quick coupling	IU 7432-30	IU 7432-30
Twist drill ø4.5, L 145, coil 50, quick coupling	IU 7445-00	IU 7445-00
Depth gauge for screws ø3.5 - 4.0, up to L 90	IS 7904-20	IS 7904-20
Screwdriver, quick coupling	IU 7835-00	IU 7835-56
Large handle, cannulated, quick coupling	IU 7706-00	IU 7706-00



◆ CAUTION:

Cortical screws are offered with different head diameters. Follow the instructions on page 15 for choosing the appropriate screw!

◆ NOTE:

If a combination of locking and non-locking screws is used, non-locking cortical screws must be inserted first.

- To insert a cortical screw (gold), place the double drill guide in the plate hole and press it down. Choose an appropriate drill and drill through both cortices. Determine the length of the screw using the depth gauge and insert a screw of appropriate length using the screwdriver.

◆ NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Check plate position using fluoroscopy and adjust if required.
- When using a cortical screw as a lag screw, start with the drill for the gliding hole and the appropriate side of the double drill guide. Drill through the near cortex or perforating the fracture line, center the other side of the drill guide in the gliding hole and drill the core hole with an appropriate diameter drill through the far cortex. Determine the screw length using the depth gauge and insert a non-locking coral screw of the appropriate length.



Insertion of cortical screws
(blue)

INSTRUMENTS

Drill guide for round hole LOQTEQ® 3.5, I-ø 2.8, blue
 Twist drill ø2.7, L 150, coil 50, quick coupling
 Twist drill ø2.7, L 150, coil 50, quick coupling, scaled
 Stop ring for depth measurement, SF
 Depth gauge for screws ø3.5-4.0, up to L 90
 Screwdriver Duo, T15, quick coupling
 Handle with quick coupling, with torque limiter, 2.0Nm

ART.-NO.

IU 8166-20
 IU 7427-15
 IU 7427-16
 IU 8166-06
 IS 7904-20
 IU 7825-56
 IU 7707-20



3.5



- Insert a drill guide ø2.7 mm (blue) into any chosen round hole and drill to the desired depth using a drill (blue/red).

◆ CAUTION:

The screwdriver duo is not intended for screwing the drill guide into the plate.

- The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.
- The stop ring facilitates reading off the calibration when attached to the drill. Push it down to the guide sleeve and remove the drill for reading the drilling depth in the gap of the ring.

◆ NOTE:

The screwdriver duo facilitates manual removal of the drill guide.

- Select a locking screw (blue) of the proper length. Loosely insert the screw using the screwdriver T15 manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface.

◆ NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Finish the screw manually using the screwdriver T15 with the torque limiting handle 2.0 Nm. Optimal locking should be achieved with an audible and tactile click of the torque limiter.

◆ CAUTION:

As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.

- Follow these instructions to insert further screws in the plate holes depending on the fracture pattern. Finally, confirm that all screw heads are flush with the plate surface. Check the result using fluoroscopy and adjust screw positioning or length as necessary.

Insertion of locking compression screws (red) without compression



INSTRUMENTS SMALL FRAGMENT

Drill guide for gliding hole LOQTEQ® 3.5, I-ø 2.8, red
Twist drill ø2.7, L 150, coil 50, quick coupling
Twist drill ø2.7, L 150, coil 50, quick coupling, scaled
Stop ring for depth measurement, SF
Depth gauge for screws ø3.5–4.0, up to L 90
Screwdriver Duo, T15, quick coupling
Handle with quick coupling, with torque limiter, 2.0Nm

ART.-NO.

IU 8166-10
IU 7427-15
IU 7427-16
IU 8166-06
IS 7904-20
IU 7825-56
IU 7707-20

INSTRUMENTS LARGE FRAGMENT

Drill guide for gliding hole LOQTEQ® 4.5, I-ø 3.9, red
Twist drill ø3.8, L 180, coil 50, quick coupling
Stop ring for depth measurement, LF
Depth gauge for screws ø4.5–6.5, up to L 100
Screwdriver Duo, T25, quick coupling
Handle with quick coupling, with torque limiter 3.5Nm

ART.-NO.

IU 8167-10
IU 7438-18
IU 8184-03
IS 7905-20
IU 7835-56
IU 7707-35



- Insert a drill guide (red) into any chosen gliding hole and drill to the desired depth using a drill (blue/red).

◆ CAUTION:

The screwdriver duo is not intended for screwing the drill guide into the plate.

- The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.
- The stop ring facilitates reading off the calibration when attached to the drill. Push it down to the guide sleeve and remove the drill for reading the drilling depth in the gap of the ring.

◆ NOTE:

The screwdriver duo facilitates manual removal of the drill guide.

- Select a locking screw (rot) of the proper length. Loosely insert the screw using the screwdriver manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface.

◆ NOTE:

Ensure proper alignment of the screwdriver and that the screwdriver tip is fully seated in the screw head.

- Finish the screw manually using the screwdriver bit with the torque limiting handle. Optimal locking should be achieved with an audible and tactile click of the torque limiter.

◆ CAUTION:

As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.

Insertion of locking compression screws (red) with compression



INSTRUMENTS SMALL FRAGMENT

Basic Insert for Load Drill Guide LOQTEQ® 3.5
 Load Drill guide LOQTEQ® 3.5, compression 1mm
 Load Drill guide LOQTEQ® 3.5, compression 2mm
 Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm

ART.-NO.

IU 8166-05
 IU 8166-01
 IU 8166-02
 IU 8166-03

INSTRUMENTS LARGE FRAGMENT

Basic insert for load drill guide LOQTEQ® 4.5
 Load Drill guide LOQTEQ® 4.5, compression 1mm
 Load Drill guide LOQTEQ® 4.5, compression 2mm
 Load drill guide LOQTEQ® 4.5, adjustable up to 2mm

ART.-NO.

IU 8167-05
 IU 8167-01
 IU 8167-02
 IU 8167-03



- LOQTEQ® gliding holes allow for fracture compression with subsequent locking fixation in one step. Load drill guides enable compressing fracture gaps of up to 2 mm.
- Screw the basic insert for load drill guides into a shaft hole near the fracture line or, if necessary, above the fracture line. Choose a load drill guide in accordance with the compression distance (1 mm or 2 mm), slide it on the basic insert and place the drill guide in the next plate hole, away from the fracture gap. Avoid pressure on the drill guide.


◆ CAUTION:

The screwdriver duo is not intended for screwing the basic insert into the plate.

- If available, use the adjustable load drill guide. The fracture gap serves as orientation in setting the compression distance (max. 2 mm) by turning the wheel of the load drill guide until an appropriate gap opens in the upper part of the instrument.

◆ NOTE:

Care should be taken when selecting the proper compression distance (1 mm or 2 mm). Avoid overcompression to ensure full locking of the screw, especially in hard bone.





INSTRUMENTS SMALL FRAGMENT

Twist drill $\varnothing 2.7$, L 150, coil 50, quick coupling
Twist drill $\varnothing 2.7$, L 150, coil 50, quick coupling, scaled
Depth gauge for screws $\varnothing 3.5$ -4.0, up to L 90
Screwdriver Duo, T15, quick coupling
Handle with quick coupling, with torque limiter, 2.0 Nm
Handle for quick coupling, large, cannulated

ART.-NO.

IU 7427-15
IU 7427-16
IS 7904-20
IU 7825-56
IU 7707-20
IU 7706-00

INSTRUMENTS LARGE FRAGMENT

Twist drill $\varnothing 3.8$, L 250, coil 50, quick coupling
Depth gauge for screws $\varnothing 4.5$ -6.5, up to L 100mm
Screwdriver Duo, T25, quick coupling
Handle with quick coupling, with torque limiter 3.5Nm
Large handle, cannulated, quick coupling

ART.-NO.

IU 7438-25
IS 7905-20
IU 7835-56
IU 7707-35
IU 7706-00

- Drill to the desired depth using a drill (blue/red) and remove the basic insert. The screw length can be read off the calibration of the drill or determined using the depth gauge, after the drill guide has been removed.

◆ NOTE:



The screwdriver duo facilitates manual removal of the basic insert.

- Select a locking screw (red) of the proper length. Loosely insert the screw using the screwdriver manually or under power with a low speed. Stop insertion when the screw head approaches the plate surface. Finish the screw manually using the screwdriver bit with the torque limiting handle. Optimal locking should be achieved with an audible and tactile click of the torque limiter.

◆ CAUTION:

As soon as the head of the screw reaches the plate hole it is compulsory to switch to the torque limiter. In cases of uncommonly hard bone, it may be necessary to finish the screw without the torque limiter to ensure the screw head is flush with the plate and the screw is locked.

- Alternatively, insert a non-locking cortical screw (gold) as a compression screw by placing the double drill guide without pressure in an off-center position in the plate hole. For insertion of the screw, follow the instructions on page 7.
- Follow these instructions to insert further screws in the plate holes depending on the fracture pattern. Finally, confirm that all screw heads are flush with the plate surface. Check the result using fluoroscopy and adjust screw positioning or length as necessary.





INSTRUMENTS

Explantation screwdriver T15, round handle
Explantation screwdriver T25, round handle

ART.-NO.


IU 7811-15
IU 7811-25

◆ NOTE:

The screwdrivers T15 (IU 7825-56) and T25 (IU 7835-56) in the set are self-retaining and should not be used for screw explantation.

- Use the corresponding explantation screwdriver for safe removal of a screw. Explantation screwdrivers are not self-retaining, penetrate further into the screw head and thus permit a higher torque when removing screws. They are not included in the set as standard and must be ordered separately.
- Place an incision on the old scar. Manually undo all screws and sequentially remove them. After manually unlocking all screws, removal may be performed using a power tool.


Disassembly



Adjustable load drill guide (IU 8166-03 / IU 8167-03)

- Remove screws (item 4.1 and 4.2) using a hexagonal screwdriver 2.5
- Unscrew the set screw (item 3)
- Pull the compression block apart (items 1 and 2)

Assembly



- Fit together the compression block (items 1 and 2)
- Insert the set screw (item 3) into the compression block, middle hole
- Insert the retaining screws (items 4.1 and 4.2) using a hexagonal screwdriver 2.5



LOQTEQ® Straight Plate 3.5

HOLES	LENGTH (mm)	ART.-NO.
4	60	PG 3555-04-2
5	73	PG 3555-05-2
6	86	PG 3555-06-2
7	99	PG 3555-07-2
8	112	PG 3555-08-2
9	125	PG 3555-09-2
10	138	PG 3555-10-2
11	151	PG 3555-11-2
12	164	PG 3555-12-2
14	190	PG 3555-14-2



LOQTEQ® Reconstruction Plate 3.5

HOLES	LENGTH (mm)	ART.-NO.
5	58	PR 3550-05-2
6	70	PR 3550-06-2
7	82	PR 3550-07-2
8	94	PR 3550-08-2
9	106	PR 3550-09-2
10	118	PR 3550-10-2
12	142	PR 3550-12-2
14	166	PR 3550-14-2
16	190	PR 3550-16-2
18	214	PR 3550-18-2



LOQTEQ® 1/3 Tubular Plate 3.5

HOLES	LENGTH (mm)	ART.-NO.
3	40	PG 3553-03-2
4	52	PG 3553-04-2
5	64	PG 3553-05-2
6	76	PG 3553-06-2
7	88	PG 3553-07-2
8	100	PG 3553-08-2
9	112	PG 3553-09-2
10	124	PG 3553-10-2
11	136	PG 3553-11-2
12	148	PG 3553-12-2

LOQTEQ® Cortical Screw 3.5, T15, self-tapping



L 12	SK 3525-12-2
L 14	SK 3525-14-2
L 16	SK 3525-16-2
L 18	SK 3525-18-2
L 20	SK 3525-20-2
L 22	SK 3525-22-2
L 24	SK 3525-24-2
L 26	SK 3525-26-2
L 28	SK 3525-28-2
L 30	SK 3525-30-2
L 32	SK 3525-32-2
L 34	SK 3525-34-2
L 36	SK 3525-36-2
L 38	SK 3525-38-2
L 40	SK 3525-40-2
L 42	SK 3525-42-2
L 44	SK 3525-44-2
L 45	SK 3525-45-2
L 50	SK 3525-50-2
L 55	SK 3525-55-2
L 60	SK 3525-60-2
L 65	SK 3525-65-2
L 70	SK 3525-70-2

LOQTEQ® Cortical Screw 3.5, small head, T15, self-tapping



L 10	SK 3526-10-2
L 12	SK 3526-12-2
L 14	SK 3526-14-2
L 16	SK 3526-16-2
L 18	SK 3526-18-2
L 20	SK 3526-20-2
L 22	SK 3526-22-2
L 24	SK 3526-24-2
L 26	SK 3526-26-2
L 28	SK 3526-28-2
L 30	SK 3526-30-2
L 32	SK 3526-32-2
L 34	SK 3526-34-2
L 36	SK 3526-36-2
L 38	SK 3526-38-2
L 40	SK 3526-40-2
L 42	SK 3526-42-2
L 45	SK 3526-45-2
L 50	SK 3526-50-2
L 55	SK 3526-55-2
L 60	SK 3526-60-2
L 65	SK 3526-65-2
L 70	SK 3526-70-2
L 75	SK 3526-75-2
L 80	SK 3526-80-2
L 85	SK 3526-85-2
L 90	SK 3526-90-2

Cortical Screw 3.5, T15, self-tapping



L 10	SK 3514-10-2
L 12	SK 3514-12-2
L 14	SK 3514-14-2
L 16	SK 3514-16-2
L 18	SK 3514-18-2
L 20	SK 3514-20-2
L 22	SK 3514-22-2
L 24	SK 3514-24-2
L 26	SK 3514-26-2
L 28	SK 3514-28-2
L 30	SK 3514-30-2
L 32	SK 3514-32-2
L 34	SK 3514-34-2
L 36	SK 3514-36-2
L 38	SK 3514-38-2
L 40	SK 3514-40-2
L 42	SK 3514-42-2
L 45	SK 3514-45-2
L 50	SK 3514-50-2
L 55	SK 3514-55-2
L 60	SK 3514-60-2
L 65	SK 3514-65-2
L 70	SK 3514-70-2
L 75	SK 3514-75-2
L 80	SK 3514-80-2
L 85	SK 3514-85-2
L 90	SK 3514-90-2

Cortical Screw 3.5, self-tapping



L 10	SK 3510-10-2
L 12	SK 3510-12-2
L 14	SK 3510-14-2
L 16	SK 3510-16-2
L 18	SK 3510-18-2
L 20	SK 3510-20-2
L 22	SK 3510-22-2
L 24	SK 3510-24-2
L 26	SK 3510-26-2
L 28	SK 3510-28-2
L 30	SK 3510-30-2
L 32	SK 3510-32-2
L 34	SK 3510-34-2
L 36	SK 3510-36-2
L 38	SK 3510-38-2
L 40	SK 3510-40-2
L 45	SK 3510-45-2
L 50	SK 3510-50-2
L 55	SK 3510-55-2
L 60	SK 3510-60-2
L 65	SK 3510-65-2
L 70	SK 3510-70-2

Cortical Screw 3.5, small head, self-tapping



L 12	SK 3512-12-2
L 14	SK 3512-14-2
L 16	SK 3512-16-2
L 18	SK 3512-18-2
L 20	SK 3512-20-2
L 22	SK 3512-22-2
L 24	SK 3512-24-2
L 26	SK 3512-26-2
L 28	SK 3512-28-2
L 30	SK 3512-30-2
L 32	SK 3512-32-2
L 34	SK 3512-34-2
L 36	SK 3512-36-2
L 38	SK 3512-38-2
L 40	SK 3512-40-2
L 45	SK 3512-45-2
L 50	SK 3512-50-2
L 55	SK 3512-55-2
L 60	SK 3512-60-2
L 65	SK 3512-65-2
L 70	SK 3512-70-2
L 75	SK 3512-75-2
L 80	SK 3512-80-2
L 85	SK 3512-85-2
L 90	SK 3512-90-2

CAUTION:

SK 3510-xx-2 **NOT** for use in round locking holes!

SK 3512-xx-2 **NOT** for use in gliding-locking holes!

Caution:

Current tray contents do no longer include these screws. Use the part numbers on the screw racks for your order or ask your local sales agent.

Cancellous Screw 4.0, small head, T15



L 10	TL 5	SP 4030-10-2
L 12	TL 5	SP 4030-12-2
L 14	TL 5	SP 4030-14-2
L 16	TL 6	SP 4030-16-2
L 18	TL 7	SP 4030-18-2
L 20	TL 8	SP 4030-20-2
L 22	TL 9	SP 4030-22-2
L 24	TL 10	SP 4030-24-2
L 26	TL 12	SP 4030-26-2
L 28	TL 14	SP 4030-28-2
L 30	TL 14	SP 4030-30-2
L 32	TL 14	SP 4030-32-2
L 34	TL 14	SP 4030-34-2
L 36	TL 14	SP 4030-36-2
L 38	TL 14	SP 4030-38-2
L 40	TL 14	SP 4030-40-2
L 42	TL 15	SP 4030-42-2
L 45	TL 15	SP 4030-45-2
L 50	TL 15	SP 4030-50-2
L 55	TL 16	SP 4030-55-2
L 60	TL 16	SP 4030-60-2
L 65	TL 16	SP 4030-65-2
L 70	TL 16	SP 4030-70-2
L 75	TL 16	SP 4030-75-2
L 80	TL 16	SP 4030-80-2
L 85	TL 16	SP 4030-85-2
L 90	TL 16	SP 4030-90-2

Cancellous Screw 4.0, kl. Kopf, T15, full thread



L 10	SP 4035-10-2
L 12	SP 4035-12-2
L 14	SP 4035-14-2
L 16	SP 4035-16-2
L 18	SP 4035-18-2
L 20	SP 4035-20-2
L 22	SP 4035-22-2
L 24	SP 4035-24-2
L 26	SP 4035-26-2
L 28	SP 4035-28-2
L 30	SP 4035-30-2
L 32	SP 4035-32-2
L 34	SP 4035-34-2
L 36	SP 4035-36-2
L 38	SP 4035-38-2
L 40	SP 4035-40-2
L 42	SP 4035-42-2
L 45	SP 4035-45-2
L 50	SP 4035-50-2
L 55	SP 4035-55-2
L 60	SP 4035-60-2
L 65	SP 4035-65-2
L 70	SP 4035-70-2
L 75	SP 4035-75-2
L 80	SP 4035-80-2
L 85	SP 4035-85-2

Washer

I-ø 4.4 mm, A-ø 8.0 mm, Titanium



SU 0448-00-2




Bending iron 1 for small fragment plates, closed
Bending iron 2 for small fragment plates, closed

IP 8405-00
IP 8405-50



Depth gauge for screws $\varnothing 3.5 - 4.0$, up to L 90

IS 7904-20



Twist drill $\varnothing 2.5$, L 110, coil 50, quick coupling
Twist drill $\varnothing 2.7$, L 150, coil 50, quick coupling
Twist drill $\varnothing 2.7$, L 150, coil 50, quick coupling, scaled
Twist drill $\varnothing 2.7$, L 220, coil 50, quick coupling
Twist drill $\varnothing 2.7$, L 220, coil 50, quick coupling, skaliert
Twist drill $\varnothing 3.5$, L 110, coil 50, quick coupling

IU 7425-00
IU 7427-15
IU 7427-16
IU 7427-22
IU 7427-23
IU 7435-00



Large handle, cannulated, quick couplin

IU 7706-00



Handle with quick coupling, with torque limiter 2.0Nm

IU 7707-20



Screwdriver Duo, T15, quick coupling
Screwdriver Duo, SW 2.5, quick coupling

IU 7825-56 ★
IU 7825-00 ●



Double drill guide $\varnothing 2.7/3.5$, with spring aided centering
Double drill guide $\varnothing 2.5/3.5$, with spring aided centering

IU 8116-60 ★
IU 8116-50 ●



Load Drill guide LOQTEQ® 3.5, compression 1mm
 Load Drill guide LOQTEQ® 3.5, compression 2mm
 Load Drill guide LOQTEQ® 3.5, adjustable up to 2mm
 Basic insert for load drill guide LOQTEQ® 3.5

IU 8166-01
 IU 8166-02
 IU 8166-03
 IU 8166-05



Stop ring for depth measurement, SF

IU 8166-06



Drill guide for gliding hole LOQTEQ® 3.5, l-ø 2.8, red

IU 8166-10




Reduction sleeve for K-wire ø1.6

IU 8166-16



Drill guide for round hole LOQTEQ® 3.5, l-ø 2.8, blue

IU 8166-20



Caddy for K-wire L 150
 K-wire with trocar point, ø1.6, L 150

IC 0006-15
 NK 0016-15





LOQTEQ® Narrow Plate 4.5

HOLES	LENGTH (mm)	ART.-NO.
4	72	PG 4555-04-2
5	90	PG 4555-05-2
6	108	PG 4555-06-2
7	126	PG 4555-07-2
8	144	PG 4555-08-2
9	162	PG 4555-09-2
10	180	PG 4555-10-2
11	198	PG 4555-11-2
12	216	PG 4555-12-2
14	252	PG 4555-14-2
16	288	PG 4555-16-2



LOQTEQ® Broad Plate 4.5

HOLES	LENGTH (mm)	ART.-NO.
6	115	PG 4556-06-2
7	133	PG 4556-07-2
8	150	PG 4556-08-2
9	168	PG 4556-09-2
10	186	PG 4556-10-2
11	204	PG 4556-11-2
12	222	PG 4556-12-2
14	257	PG 4556-14-2
16	293	PG 4556-16-2
18	328	PG 4556-18-2

**LOQTEQ® Kortikalisschraube 4.5,
T25, self-tapping**



L 14	SK 4525-14-2
L 16	SK 4525-16-2
L 18	SK 4525-18-2
L 20	SK 4525-20-2
L 22	SK 4525-22-2
L 24	SK 4525-24-2
L 26	SK 4525-26-2
L 28	SK 4525-28-2
L 30	SK 4525-30-2
L 32	SK 4525-32-2
L 34	SK 4525-34-2
L 36	SK 4525-36-2
L 38	SK 4525-38-2
L 40	SK 4525-40-2
L 42	SK 4525-42-2
L 45	SK 4525-45-2
L 50	SK 4525-50-2
L 55	SK 4525-55-2
L 60	SK 4525-60-2
L 65	SK 4525-65-2
L 70	SK 4525-70-2
L 75	SK 4525-75-2
L 80	SK 4525-80-2
L 85	SK 4525-85-2
L 90	SK 4525-90-2

**Kortikalisschraube 4.5,
small head, self-tapping***



L 14	SK 4512-14-2
L 16	SK 4512-16-2
L 18	SK 4512-18-2
L 20	SK 4512-20-2
L 22	SK 4512-22-2
L 24	SK 4512-24-2
L 26	SK 4512-26-2
L 28	SK 4512-28-2
L 30	SK 4512-30-2
L 32	SK 4512-32-2
L 34	SK 4512-34-2
L 36	SK 4512-36-2
L 38	SK 4512-38-2
L 40	SK 4512-40-2
L 42	SK 4512-42-2
L 45	SK 4512-45-2
L 50	SK 4512-50-2
L 55	SK 4512-55-2
L 60	SK 4512-60-2
L 65	SK 4512-65-2
L 70	SK 4512-70-2
L 75	SK 4512-75-2
L 80	SK 4512-80-2
L 85	SK 4512-85-2
L 90	SK 4512-90-2

Washer

I-ø 5.2 mm, A-ø 10.0 mm, Titanium



SU 0510-00-2

**Cortical Screw 4.5,
T25, self-tapping**



L 20	SK 4514-20-2
L 22	SK 4514-22-2
L 24	SK 4514-24-2
L 26	SK 4514-26-2
L 28	SK 4514-28-2
L 30	SK 4514-30-2
L 32	SK 4514-32-2
L 34	SK 4514-34-2
L 36	SK 4514-36-2
L 38	SK 4514-38-2
L 40	SK 4514-40-2
L 42	SK 4514-42-2
L 45	SK 4514-45-2
L 50	SK 4514-50-2
L 55	SK 4514-55-2
L 60	SK 4514-60-2
L 65	SK 4514-65-2
L 70	SK 4514-70-2
L 75	SK 4514-75-2
L 80	SK 4514-80-2
L 85	SK 4514-85-2
L 90	SK 4514-90-2

**Cortical Screw 4.5,
self-tapping***



L 14	SK 4510-14-2
L 16	SK 4510-16-2
L 18	SK 4510-18-2
L 20	SK 4510-20-2
L 22	SK 4510-22-2
L 24	SK 4510-24-2
L 26	SK 4510-26-2
L 28	SK 4510-28-2
L 30	SK 4510-30-2
L 32	SK 4510-32-2
L 34	SK 4510-34-2
L 36	SK 4510-36-2
L 38	SK 4510-38-2
L 40	SK 4510-40-2
L 42	SK 4510-42-2
L 45	SK 4510-45-2
L 50	SK 4510-50-2
L 55	SK 4510-55-2
L 60	SK 4510-60-2
L 65	SK 4510-65-2
L 70	SK 4510-70-2

Caution:

SK 4510-xx-2 **NOT** for use in round locking holes!

*** Caution:**

Current tray contents do no longer include these screws. Use the part numbers on the screw racks for your order or ask your local sales agent.



Depth gauge for screws $\varnothing 4.5 - 6.5$, up to L 100

IS 7905-20



Twist drill $\varnothing 3.2$, L 195, coil 50, quick coupling

IU 7432-30

Twist drill $\varnothing 3.8$, L 180, coil 50, quick coupling

IU 7438-18

Twist drill $\varnothing 3.8$, L 250, coil 50, quick coupling

IU 7438-25

Twist drill $\varnothing 4.5$, L 145, coil 50, quick coupling

IU 7445-00



Large handle, cannulated, quick coupling

IU 7706-00



Large handle, cannulated, quick coupling

IU 7707-35



Screwdriver Duo, T25, quick coupling

IU 7835-56 ★

Screwdriver, hexagonal, $\varnothing 3.5$ for quick coupling

IU 7835-00 ★



Double drill guide $\varnothing 3.2/4.5$, with spring aided centering

IU 8117-50



Load drill guide LOQTEQ® 4.5, compression 1mm
Load drill guide LOQTEQ® 4.5, compression 2mm
Load drill guide LOQTEQ® 4.5, adjustable up to 2mm
Basic insert for load drill guide LOQTEQ® 4.5

IU 8167-01
IU 8167-02
IU 8167-03
IU 8167-05



Bohrbuchse für Gleitloch LOQTEQ® 4.5, I-ø 3.9, rot

IU 8167-10




Reduction sleeve for K-wire ø2.0

IU 8167-15



Stop ring for depth measurement, LF

IU 8184-03




Caddy for K-wire L 250
K-wire with trocar point, ø2.0, L 250

IC 0006-25*
NK 0020-25

* Not included in the set, must be ordered separately

Preoperative




Postoperative




Clinical case and CT images with the kind permission of Asklepios Clinic Weißenfels, Germany

Ulna shortening osteotomy

Preoperative



Postoperative



Clinical case and CT images with the kind permission of University Hospital Gießen and Marburg, Germany

Subject to technical modifications,
errors and misprints.

© aap Implantate AG
WP 40P001 EN / 2011-2

aap Implantate AG
Lorenzweg 5 • 12099 Berlin
Germany

Phone +49 30 75019-0
Fax +49 30 75019-111

customer.service@aap.de
www.aap.de

LOCTEC®



aap Implantate AG
Lorenzweg 5 • 12099 Berlin
Germany

Phone +49 30 75019-0
Fax +49 30 75019-111

customer.service@aap.de
www.aap.de



(01)04042409437182(10)2011-2
WP 40P001 EN / 2011-2