



LOQTEQ® Ellenbogenplatten 2.7/3.5 • Elbow Plates 2.7/3.5

IC 6935-00 Implantateset • Set of Implants



ARTIKEL	ARTICLE	ANZAHL QUANTITY	ART.-NR. ART.-NO.
Leersieb Implantate, LOQTEQ® Ellenbogenplatten 2.7/3.5	Tray for instruments, LOQTEQ® Elbow, empty	1	IC 6935-01
Deckel für Leersiebe, groß	Lid for trays, large	1	IC 2008-00
1 Messinstrument für Schrauben ø2.7, bis L 70	Depth gauge for screws ø2.7, up to L 70	1	IS 7903-20
2 Bohrer Schnellkupplung ø2.0, L 180, Wendel 50	Twist drill ø2.0, L 180, coil 50, quick coupling	2	IU 7420-18
3 Handgriff rund für Schnellkupplung mit Drehmomentbegrenzer 1.5 Nm	Handle round with quick coupling, with torque limiter 1.5 Nm	1	IU 7707-00
4 Schraubendrehereinsatz Duo, T8, Schnellkupplung	Screwdriver Duo, T8, quick coupling	2	IU 7815-56
5 Doppelbohrbuchse, Bohrer ø2.0/2.5	Double drill guide ø2.0/2.5	1	IU 8125-00
6 Bohrbuchse LOQTEQ® 2.7, Skala bis L 75, Bohrer ø2.0, hellblau	Drill guide LOQTEQ® 2.7, scale to L 75, drill ø2.0, light blue	2	IU 8169-20
7 Zielblock LOQTEQ® Olekranonplatte, R	Aiming device LOQTEQ® Olecranon Plate, R	1	IU 8178-01
8 Zielblock LOQTEQ® Olekranonplatte, L	Aiming device LOQTEQ® Olecranon Plate, L	1	IU 8178-02
9 Zielbügel LOQTEQ® Distale Mediale Humerusplatte	Aiming arm LOQTEQ® Distal Medial Humerus Plate	1	IU 8179-00
LOQTEQ® Distale Mediale Humerusplatte, 2-Loch, L 85, R, Titan	LOQTEQ® Distal Medial Humerus Plate, 2 holes, L 85, R, Titanium	1	PH 3521-02-2
10 LOQTEQ® Distale Mediale Humerusplatte, 3-Loch, L 98, R, Titan	LOQTEQ® Distal Medial Humerus Plate, 3 holes, L 98, R, Titanium	1	PH 3521-03-2
LOQTEQ® Distale Mediale Humerusplatte, 5-Loch, L 124, R, Titan	LOQTEQ® Distal Medial Humerus Plate, 5 holes, L 124, R, Titanium	1	PH 3521-05-2
LOQTEQ® Distale Mediale Humerusplatte, 2-Loch, L 85, L, Titan	LOQTEQ® Distal Medial Humerus Plate, 2 holes, L 85, L, Titanium	1	PH 3522-02-2
11 LOQTEQ® Distale Mediale Humerusplatte, 3-Loch, L 98, L, Titan	LOQTEQ® Distal Medial Humerus Plate, 3 holes, L 98, L, Titanium	1	PH 3522-03-2
LOQTEQ® Distale Mediale Humerusplatte, 5-Loch, L 124, L, Titan	LOQTEQ® Distal Medial Humerus Plate, 5 holes, L 124, L, Titanium	1	PH 3522-05-2
LOQTEQ® Distale Dorsolat. Humerusplatte, 2-Loch, L 89, R, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 2 holes, L 89, R, Titanium	1	PH 3531-02-2
12 LOQTEQ® Distale Dorsolat. Humerusplatte, 3-Loch, L 102, R, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 3 holes, L 102, R, Titanium	1	PH 3531-03-2
LOQTEQ® Distale Dorsolat. Humerusplatte, 5-Loch, L 128, R, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 5 holes, L 128, R, Titanium	1	PH 3531-05-2
LOQTEQ® Distale Dorsolat. Humerusplatte, 2-Loch, L 89, L, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 2 holes, L 89, L, Titanium	1	PH 3532-02-2
13 LOQTEQ® Distale Dorsolat. Humerusplatte, 3-Loch, L 102, L, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 3 holes, L 102, L, Titanium	1	PH 3532-03-2
LOQTEQ® Distale Dorsolat. Humerusplatte, 5-Loch, L 128, L, Titan	LOQTEQ® Distal Dorsolat. Humerus Plate, 5 holes, L 128, L, Titanium	1	PH 3532-05-2
LOQTEQ® Olekranonplatte, 2-Loch, L 86, R, Titan	LOQTEQ® Olecranon Plate, 2 holes, L 86, R, Titanium	1	PU 3531-02-2
14 LOQTEQ® Olekranonplatte, 4-Loch, L 112, R, Titan	LOQTEQ® Olecranon Plate, 4 holes, L 112, R, Titanium	1	PU 3531-04-2
LOQTEQ® Olekranonplatte, 2-Loch, L 86, L, Titan	LOQTEQ® Olecranon Plate, 2 holes, L 86, L, Titanium	1	PU 3532-02-2
15 LOQTEQ® Olekranonplatte, 4-Loch, L 112, L, Titan	LOQTEQ® Olecranon Plate, 4 holes, L 112, L, Titanium	1	PU 3532-04-2
Schraubenbank leer, Erweiterungsset LOQTEQ® Ellenbogen	Screw rack for extension set LOQTEQ® Elbow, empty	1	IC 6935-02
16 Kortikalisschraube 2.5, kl. Kopf T8, selbstschn., L 14-70, Titan	Cortical Screw 2.5, small head T8, self-tapp., L 14-70, Titanium	4	SK 2512-xx-2
17 LOQTEQ® Kortikalisschr. 2.7, kl. Kopf T8, selbstschn., L 14-70, Titan	LOQTEQ® Cortical Screw 2.7, small head T8, self-tapp., L 14-70, Titanium	8	SK 2726-xx-2