



LOQTEQ® Osteotomie Set, HTO / DFO • Osteotomy Set HTO / DFO

IC 4500-30



ARTIKEL	ARTICLE	ANZAHL QUANTITY	ART.-NR. ART.-NO.
Leersieb für LOQTEQ® Osteotomie System	Tray for LOQTEQ® Osteotomy System, empty	1	IC 4500-31
Deckel für Leersiebe, groß	Lid for trays, large	1	IC 2008-00
1 Messinstrument für Schrauben ø4.5-6.5, bis L 100mm	Depth gauge for screws ø4.5-6.5, up to L 100mm	1	IS 7905-20
2 Knochenspreizzange, L 280	Bone spread forceps, L 280	0	IU 2503-00
3 Lambotte-Meißel, 15mm	Lambotte chisel, 15mm	1	IU 3000-15
4 Lambotte-Meißel, 20mm	Lambotte chisel, 20mm	2	IU 3000-20
5 Bohrer Schnellkupplung ø3.2, L 195, Wendel 50	Twist drill ø3.2, L 195, coil 50, quick coupling	1	IU 7432-30
6 Bohrer Schnellkupplung ø3.8, L 200, Wendel 50	Twist drill ø3.8, L 200, coil 50, quick coupling	1	IU 7438-20
7 Handgriff für Schnellkupplung groß, durchbohrt	Large handle, cannulated, quick coupling	1	IU 7706-00
8 Handgriff für Schnellkupplung m. Drehmomentbegrenzer 3.5 Nm	Handle with quick coupling, with torque limiter 3.5 Nm	1	IU 7707-35
9 Schraubendrehereinsatz Duo, T25, Schnellkupplung	Screwdriver Duo, T25, quick coupling	2	IU 7835-56
10 Messstab für K-Draht ø2.5, L 200	Measuring device for K-wire 2.5, L 200	1	IU 7925-20
11 Keilmesslehre	Wedge gauge	1	IU 7960-00
12 Winkellehre für Closed-Wedge-Osteotomie	Angle gauge for closed wedge osteotomy	0	IU 7970-00
13 Abstandshalter Winkelstabil LOQTEQ® 4.5	Angle stable locking spacer LOQTEQ® 4.5	4	IU 7972-00
13 Abstandshalter Winkelstabil LOQTEQ® 4.5	Angle stable locking spacer LOQTEQ® 4.5	0	IU 7972-00
14 Externe Ausrichtung	External alignment device	0	IU 7973-00
15 Spannbohrbuchse LOQTEQ® 4.5, variabel bis 2mm	Load drill guide LOQTEQ® 4.5, adjustable up to 2mm	0	IU 8167-03
16 Reduzierhülse für K-Draht ø2.0	Reduction sleeve for K-wire ø2.0	2	IU 8167-15
17 Bohrbuchse für Rundloch LOQTEQ® 4.5, I-ø 3.9, blau	Drill guide for round hole LOQTEQ® 4.5, I-ø3.9,3 blue	3	IU 8167-20
18 Grundhalterung für Spannbohrbuchse, LOQTEQ® 4.5, Rundloch	Basic insert for load drill guide, LOQTEQ® 4.5, round hole	0	IU 8167-45
19 Befestigungsschraube Zielblock LOQTEQ® GFI T25	Fixing screw aiming device LOQTEQ® LFI T25	3	IU 8176-04
20 Zielblock LOQTEQ® HTO-Platte	Aiming device LOQTEQ® HTO plate	1	IU 8184-01
21 Anschlagring für Tiefenmessung, GF	Stop ring for depth measurement, LF	1	IU 8184-03
22 Zielblock LOQTEQ® DFO-Platte, R	Aiming device LOQTEQ® DFO plate, R	1	IU 8185-01
23 Zielblock LOQTEQ® DFO-Platte, L	Aiming device LOQTEQ® DFO plate, L	1	IU 8185-02
24 Parallelbohrlehre für K-Draht, ø2.5, schmal, PEEK	Parallel drill guide for K-wire, ø2.5, narrow, PEEK	1	IU 8188-25
25 LOQTEQ® Führungshülse 4.5, blau	LOQTEQ® Screw guide sleeve 4.5, blue	1	IU 8210-45
26 K-Draht mit Trokarspitze, ø2.0, L 250	K-wire with trocar point, ø2.0, L 250	10	NK 0020-25
27 K-Draht mit Trokarspitze, ø2.5, L 200	K-wire with trocar point, ø2.5, L 200	10	NK 0025-20
28 Behälter für K-Drähte L 200	Caddy for K-wire L 200	1	IC 0006-20
29 Behälter für K-Drähte L 250	Caddy for K-wire L 250	1	IC 0006-25
30 LOQTEQ® High Tibia Osteotomieplatte 4.5, Titan	LOQTEQ® High Tibia Osteotomy Plate 4.5, Titanium	1	PO 4560-01-2
31 LOQTEQ® Distale Femur Osteotomieplatte 4.5, R, Titan	LOQTEQ® Distal Femur Osteotomy Plate 4.5, R, Titanium	1	PO 4561-01-2
32 LOQTEQ® Distale Femur Osteotomieplatte 4.5, L, Titan	LOQTEQ® Distal Femur Osteotomy Plate 4.5, L, Titanium	1	PO 4562-01-2
33 Schraubenbank leer, Erweiterungsset LOQTEQ® Osteotomie	Screw rack for extension set LOQTEQ® Osteotomy, empty	1	IC 4500-04
34 Kortikalisschraube 4.5, T25, selbstschn. L 30-50, Titan	Cortical Sreew 4.5, T25, self-tapping, L 30-50, Titanium	2	SK 4514-xx-2
35 LOQTEQ® Kortikalisschraube 4.5, T25, selbstschn. L 18-24, Titan	LOQTEQ® Cortical Screw 4.5, T25, self-tapping, L 18-24, Titanium	2	SK 4525-xx-2
36 LOQTEQ® Kortikalisschr. 4.5, kl. Kopf, T25, selbstschn. L 18-45, Titan	LOQTEQ® Cortical Screw 4.5, small head, T25, self-tapping, L 18-45, Titanium	3	SK 4526-xx-2
37 LOQTEQ® Kortikalisschr. 4.5, kl. Kopf, T25, selbstschn. L 50-85, Titan	LOQTEQ® Cortical Screw 4.5, small head, T25, self-tapping, L 50-85, Titanium	4	SK 4526-xx-2